
RM310 19 INCH RACK KIT – QPX Series

The QPX Series rack mount kit is designed to accommodate one instrument centrally in a standard 19" wide rack. Two vent panels are provided for fitting to the tray either side of the instrument.

It is recommended that side supports, e.g. slides or angle supports, are used with the tray; the tray should be drilled and assembled with the side supports before the instrument is fitted to it.

Turn the instrument upside down so that it is laying on its top face and remove the 4 screws which hold the feet; the screws retaining the tilting front feet are accessed by first prising out the self-adhesive rubber inserts in the moulded feet.

Having first fitted the vent panels to the sides of the tray using 4 of the M4 countersunk screws provided, lay the tray over the instrument with the front edge aligned with the moulded bezel of the instrument.

The screw-holes for the rear feet and 2 screw-holes near the front should line up with the appropriate holes in the tray; secure the instrument to the tray with 4 of the M4 countersunk screws provided.

Turn over the tray and fit to the rack enclosure.

Thurlby Thandar Instruments Ltd., Glebe Road, Huntingdon, Cambs PE29 7DR

Instruction Leaflet 48511-1380 Issue 2
